



## Refill kit for 32313. CO2 calibration starter kit

**Termékszám:** 32314

**ELAVULT CIKK**

Ez a termék már nem elérhető.

Reagents to prepare a reference solution to calibrate CO2 analyzers. High carbonation 3 to 8 g/kg. Refill kit for 32313.

### Traceable and Trustworthy Results

Effective calibration is critical to the performance of any analytical instrument and the ORBISPHERE 3231X is designed to give a traceable and trustworthy calibration standard with defined carbon dioxide content to +/- 0.5 %. Each reagent kit contains a certificate of analysis.

### Quick and Easy to Use

---

### Műszaki adatok

|                      |                                          |
|----------------------|------------------------------------------|
| Méréstartomány:      | 3 - 8 g/kg                               |
| Műszer:              | Orbisphere 3658 Analyzer Portable CO2    |
| Paraméter:           | CO2                                      |
| Tárolási feltételek: | Store in a dry place at room temperature |