



Beaker, Griffin, Low Form, Glass, 600 mL

Termékszám:	50052
HUF Ár (nettó):	Elérhetőségünk
Szállítás 2 héten belül	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Műszaki adatok

Anyag:	Glass, Griffin, Low Form
Anyaga:	Glass
Egység:	1 each
Kapacitás:	600 mL