



QL filters, pore size 8 - 12 μ m, diameter 70 mm, circular, 100/ pk, qualitative

Termékszám: 50653

ELAVULT CIKK
Ez a termék már nem elérhető.

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Műszaki adatok

Anyag:	Cellulose
Átmérő:	70 mm
Térfogat/ csomagméret:	100 /pk