



## QL filters, pore size 8 - 12 $\mu$ m, diameter 90 mm, circular, 100/pk, qualitative

|                              |                |
|------------------------------|----------------|
| Termékszám:                  | 50655          |
| HUF Ár (nettó):              | Elérhetőségünk |
| Szállítási dátum nem látható |                |

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Műszaki adatok

|                        |           |
|------------------------|-----------|
| Anyag:                 | Cellulose |
| Átmérő:                | 90 mm     |
| Térfogat/ csomagméret: | 100 /pk   |