



**Be Right™**



## **QL filters, pore size 8 - 12 µm, diameter 150 mm, circular, 100/pk, qualitative**

**Termékszám:**

**50658**

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### **Műszaki adatok**

Anyag: Cellulose

Átmérő: 150 mm

Térfogat/ csomagméret: 100 /pk