



Be Right™



QL filters, pore size 8 - 12 µm, diameter 185 mm, circular, 100/pk, qualitative

Termékszám:

50659

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Műszaki adatok

Anyag: Cellulose

Átmérő: 185 mm

Térfogat/ csomagméret: 100 /pk