



KNF300 pump 6 Month Service Kit

Termékszám:	19-KIT-110
HUF Ár (nettó):	Elérhetőségünk
Elérhető	

For BioTector Series 4 & B7000 analysers supplied before July 2014 – service kit to be used every 6 months for the routine maintenance of the KNF 300 circulation pump

Műszaki adatok

A doboz tartalma:	Includes: PTFE Diaphragm and NF300 O-Ring
-------------------	---

A doboz tartalma

Includes: PTFE Diaphragm and NF300 O-Ring