



B4M TOC Analyser 6 Month Service Kit

Termékszám:

19-KIT-125

Veszélyes

HUF Ár (nettó):

Elérhetőségünk

Elérhető



For BioTector B7000 MR analysers supplied from July 2014 to August 2015 with

analyser serial number format BAC-4xxx-xxxx

Műszaki adatok

A doboz tartalma:

Includes: Omron Relays (3), EMPP Tubing, Y Tube fitting (2), SR25 pump tube rail (2), PTFE Diaphragm, CO₂ Filter, Ozone Destructor catalyst, Filter Mats for vent (2), Lens Tissue for CO₂ Analyser, FKM/FPM o-ring (2)

A doboz tartalma

Includes: Omron Relays (3), EMPP Tubing, Y Tube fitting (2), SR25 pump tube rail (2), PTFE Diaphragm, CO₂ Filter, Ozone Destructor catalyst, Filter Mats for vent (2), Lens Tissue for CO₂ Analyser, FKM/FPM o-ring (2)