



**Be Right™**



## **De-bubbler including tubing and pressure relief valves for ScrubberMonitor cabinets**

**Termékszám:** **APPS10011888**

HUF Ár (nettó): Elérhetőségünk

Szállítási dátum nem látható

Complete de-bubbler kit for the ScrubberMonitor EGCS washwater monitoring cabinets. Made with marine-grade PVC built to withstand harsh environments, this system is designed to take air out of the incoming sample to provide consistent sample flow to the ScrubberMonitor. Includes all of the piping, pressure regulators, and flexible tubing needed to connect to the cabinet.