



Reagent Kit for EZ4041 Total Hardness Analyser (acid, buffer, colour and EDTA solution), 10-250 mg/L, 10-25% range

Termékszám:

APPC4041KTOL-EU **Veszélyes**

HUF Ár (nettó):

Elérhetőségünk

Szállítási dátum nem látható



Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Műszaki adatok

A doboz tartalma:

APPC4041-01-EU Acid Solution

APPC4041-02-EU Buffer Solution

APPC4041-03 Colour Solution

APPC4041-04 EDTA Solution (0.005 M)

Paraméter:

Total Hardness

A doboz tartalma

APPC4041-01-EU Acid SolutionAPPC4041-02-EU Buffer SolutionAPPC4041-03 Colour SolutionAPPC4041-04 EDTA Solution (0.005 M)